

nephrolink

Nephrology Handheld Billing Application



nephrolink handheld application makes it possible for doctors to manage and maintain patients data easily with more accuracy and agility

The nephrolink handheld application makes it possible for the doctors to manage and maintain patients' data easily with more accuracy and agility. The user friendly interface and synchronization process makes it the only portable software package that can help doctors to collect and manage patient data just by few clicks to maintain records, billing and other purposes.

The handheld application can specifically be used by doctors to:

- Collect, manage and edit patient data
- Record patients hospital visits
- Record patients dialysis visits
- Record patients diagnosis
- Sync data with online server

Practices can use this application to:

- Manage doctors associated with practice
- Manage patients history
- Manage Doctor's visits
- Import/Export visits data
- Calculate billing information

The data collection is primarily used for record and billing purposes and all collected data is being maintained on an online server. Doctors can sync their handheld device data anytime with the server; using wireless internet connection (WiFi or WiMax), GPRS connection or PC internet connection via USB. Web based application provides comprehensive interface for practices to manage all of their data. The doctor's visits data can be used for billing purposes.

Nephrolink provides mobile productivity

In hospitals, the doctors use paper pads to record patients' data on daily basis. This technique has so many vulnerabilities like:

- Costly to produce
- Prone to errors
- Usually outdated
- Record management issues
- Difficult to sum up visits by end of the month
- Data archiving issues

Nephrolink is an affordable alternative which makes it easy to collect multitudes of patients' information and visits with accuracy while on the move.




Nephrolink Features

Nephrolink provides comprehensive features to facilitate doctors in their daily tasks. Its PDA application provides following functionalities:

- Patients record management; add, edit and view patients
- Data sorting and filtering by hospitals, facilities, shifts, doctors and dates
- Hospital and Dialysis visits creation
- Auto daily and monthly visits creation
- Patients diagnosis recording and history
- Patients sharing; doctors can share patients with other doctors to perform necessary operations in case of any non-availability.
- Automatic Application update; this feature enables the application to update itself, by triggering this feature, application automatically searches for any latest updates.
- Server Sync & Data Submitting; this feature is responsible for the synchronization of data to the online web server.

Its web application enables Practices to manage their all data online without any hassle. It provides following features:

- Basic reference data management; CPT Codes, ICD9 Codes, Facilities, Hospitals, Shifts etc.
- Add/View Patients
- Import/Export data in .xls file format
- View hospital/dialysis visits
- Create/view practice users
- Define practice user access rights
- Create/view doctors



Nephrolink is not only responsible for substantial cost savings, but its major objective is to develop a sense of trust in its users regarding the secure information handling. Also, Nephrolink solution is so easy to use, even non-technical types feel comfortable with it.

Nephrolink System Requirements

PDA Application

- Windows Mobile or Windows CE 4.0/5.0 enable Pocket PC
- Supported CPU chips: Intel(R) PXA272, ARM based
- RAM: 64MB minimum
- Internet Connection; WiFi/WiMax or GPRS or PC internet connection via USB
- Microsoft ActiveSync 4.5 or higher

Web Application

- Windows Vista/Windows XP Professional or Home Edition/Windows 2003/Windows 2000
- Microsoft Internet Explorer 6.0 or above/Mozilla FireFox 2.0.0.3 or above/Opera 9.0 or above